

Urine Analyzer



- ▶ Multi-Items in 1 Test
- ▶ i-Cloud System
- ▶ Professional/Self Test

- ▶ Single Channel
- ▶ 1 Minute Rapid Test
- ▶ Intelligent & Portability

Urine Analyzer


SPEEDRY UR APPLICATION SCENARIO



**Clinical
Privacy Hospital**



**Chronic Disease
Agency**



Reproductive Hospital



Family Test

SPEEDRY UR PRODUCT PERFORMANCE PARAMETERS

Testing Principle	Light Reflection Colorimetry
Measured Item	11 /14 in one test
Measurement Pattern	Single Channel
Stored Value	1000 Results
Battery	1200mA
LED	4.3-inch color LCD with touch function
Adapter	5V , 2A;100V-240V ± 10%, 50/60Hz
Work Environment	Temperature: 5°C- 35°C; RH: ≤80%; Atmospheric Pressure: 70kPa~106kPa
Storage Environment	Temperature: -20°C- 55°C; RH: ≤85%; Atmospheric Pressure: 70kPa~106kPa
Body Size	80mmx144mmx32mm
Service Life	Five Years

SPEEDRY UR MAIN TEST ITEM

Test Item	Reference Range
PH (PH)	PH4.5-8.0
Protein(PRO)	Negative
Bilirubin (BIL)	Negative
Ketone Body (KET)	Negative
Occult Blood or Red Blood Cells (BLD)	Negative
Vitamin C(VitC)	Negative
Urinary Calcium(Ca)	Negative

Test Item	Reference Range
Specific Gravity (SG)	1.015-1.025
Glucose (GLU)	Negative
Urobilinogen (URO)	Negative or Weakly Positive
Nitrite (NIT)	Negative
Leukocyte (LEU)	Negative
Urine Creatinine(Cr)	Negative
Urine Microalbumin(UMA)	Negative

Analog Test and Measurement Philippines Inc.
 Contact Numbers: 0917-524-0067/0917-848-0857
 (02) 8881-9669
 sales@ehealthcareinc.com
<https://ehealthcareinc.com>
<https://www.facebook.com/eHealthcareInc>

